

CEREC

A QUICK PROCEDURE
WITH A LONG TRADITION

The CEREC procedure was developed at Zurich University and has been successfully used in dental practices since 1985 - for 30 years.

Worldwide more than 45,000 practices in over 50 countries utilize the CEREC system.

Since 1985 over 8 million patients have been treated using the CEREC procedure. More than 20 million ceramic restorations have been created and placed.


Long-term clinical studies demonstrate - in terms of durability - CEREC restorations are at least equivalent to gold restorations.

Depending on the individual case, CEREC restorations can achieve even better results than gold ones.

CEREC

BENEFITS AT A GLANCE

- ✓ One Visit
- ✓ No Impressions
- ✓ No Temporaries
- ✓ Metal-Free
- ✓ Highly Aesthetic
- ✓ Long Lasting



MODERN CAD/CAM TECHNOLOGY
MAKES IT POSSIBLE.

Subject to technical changes and errors in the text due to client distribution. A feature of the Cerec Academy www.cerectoctors.com | V.1 2014

ASK OUR OFFICE HOW
CEREC
CAN BENEFIT YOU!

A BEAUTIFUL SMILE WITH
CEREC



CERAMIC RESTORATIONS WITH CEREC.



CEREC RESTORATIONS CAN BE CREATED AND PLACED DURING A SINGLE APPOINTMENT.

CEREC is a sophisticated CAD/CAM system for the production of all-ceramic inlays, onlays, partial crowns, veneers and crowns for anterior and posterior teeth.

THE CEREC PROCESS 5 SIMPLE STEPS TO YOUR NEW TOOTH




STEP 1

Before treatment begins we inform you about the benefits of CEREC.



STEP 2

After removing the decayed tissue or the defective filling we create the 3D image of the tooth with the CEREC camera.




STEP 3

The restoration is designed in office on the CEREC while you wait.



STEP 4

The CEREC milling unit makes the restoration out of a ceramic block.



STEP 5

The milled restoration is bonded directly to your tooth.

THE ADVANTAGES OF CEREC

NO IMPRESSIONS

The CEREC 3D measuring camera scans the prepared tooth in just a few seconds. This eliminates the need for unpleasant impressions.

NO TEMPORARIES

CEREC restorations are milled out of a solid ceramic block in just a few minutes. They are ready to be placed immediately. This eliminates the need for temporary fillings and crowns. And you don't have to wait for a days on end until the restoration is finished.

NO LONG WAIT PERIODS

Because CEREC restorations can be placed immediately, you need only one appointment. And you benefit from tooth restorations which have demonstrated their durability time and time again.

CERAMIC BLOCK

MILLED TOOTH

BLOCK EXCESS (REMOVED)

BLOCK BASE - HOLDS BLOCK (REMOVED AFTER MILLING)